

Digital Storage Oscilloscope [DSO] Presentation based on the TDS210/TDS1002 Series user manual

Jerin Joe

Wadhvani Electronics Laboratory
Indian Institute of Technology Bombay, Mumbai, India

July 7, 2016



Oscilloscope

Oscilloscope is an electronic measuring instrument which is used to display and analyze various waveforms.

There are two types of oscilloscopes based on the mode of operation.

- Analog Oscilloscope
- Digital Storage Oscilloscope[DSO]

Cathode Ray Oscilloscope




Figure: Front panel of CRO

A CRO uses analog circuits to sense and display information on a screen. Typically, it does not have the ability to store information about captured waveforms.

Working of CRO




Figure: Working Principle of CRO

- Depicting the cathode ray tube which emits electrons from the (filament) cathode.
- The emitted electrons is then accelerated towards the screen and strikes on the fluorescent screen, due to its fluorescent properties it glows when electron strikes on it.
- The amplitude and frequency of the signal displayed on the screen depends on the voltages on vertical and horizontal deflection plates respectively.

Digital Storage Oscilloscope




Figure: Digital Storage Oscilloscope

A DSO typically digitizes and stores captured signals in an internal memory. The stored signals can then be analyzed and displayed on a screen. This offers great flexibility in analysis of the captured waveforms.

Working of DSO




Figure: Block diagram of DSO

- DSO's first (input) stage is a vertical amplifier. Vertical controls allow you to adjust the amplitude and position range at this stage.
- ADC converts the signal's voltage at these points into digital values called sample points.
- The DSO's signal path includes a microprocessor through which the measured signal passes on its way to the display.

Front Panel



Figure: The front panel of TDS210/1002B

Some other DSOs available in the lab are TDS1002, TDS1002B and GDS-1072A-U- GW INSTRUK

Checking the working of DSO

To verify the instrument is working properly.

- Take the signal probe (Eg: TDS2200). This probe has two attenuation settings, 1X and 10X. Set the attenuation button at the tip of the probe to 1X.
- Now, attach the probe tip to the PROBE COMP 5V connector and the probe reference lead to the PROBE COMP ground connector.
- Plug the probe into Channel 1 on the oscilloscope and push the AU-TOSET2 button. On the display you should see a square wave of approximately 5V peak-to-peak at 1KHz,

Triggering

The trigger determines when the oscilloscope starts to acquire data and display a waveform. When a trigger is set up properly, it can convert unstable displays or blank screens into meaningful waveforms




Figure: Untriggered Waveform





Figure: Triggered Waveform



- **Level and Hold Off** : Level sets the amplitude level and hold off sets the time before another trigger.
- **trigger Menu** : Displays the trigger menu.
- **Set level to 50%** Sets level to 50% of signal level.
- **Force trigger**: Starts an aquisition regardless of an adequate data.
- **Trigger View** : Displays trigger waveform in place of channel waveform when trigger knob is held.

Figure: Triggering controls

Capturing A signal using a DSO

Assemble the circuit as shown in the figure




Figure: Test circuit

Connect the output of our circuit (node 2) to Channel 1 of the DSO.
Are you able to see the signal properly on DSO?

To measure signal frequency, period peak-to-peak amplitude etc. of the output waveform follow the steps given below

- Push the *AUTOSET* button automatically sets the trigger levels and displays the waveforms
- Push the *MEASURE* button to go to the Measure menu
- Push the top menu box button to select Source
- Select *CH1* for the first three measurements.
- Push the top menu box button to select *Type*
- Push the first CH1 menu box button to select *Freq*
- Push the second CH1 menu box button to select *Period*
- Push the third CH1 menu box button to select *Pk-Pk*

Use cursors to take measurements

- Push the *CURSOR* button to go to the Cursor menu
- Push the top menu box button and select *Type* as *Voltage*¹
- Set the *Source* option to *CH1*
- Vary the position of the two horizontal cursor lines using the *CURSOR 1* and *CURSOR 2* knobs
- Observe the voltage values corresponding to cursors under *Cursor 1* and *Cursor 2* the difference between them under *Delta*
- Push the top menu box button and select *Type* as *Time*
- Vary the cursors using the previous knobs and observe the time values

¹Two horizontal lines will appear if you set *Type* as *Voltage* and two vertical lines will show up if you set it as *Time*

There are two channels in the DSO which can be controlled independently

- Take another probe and plug it to the *Channel 2* on the oscilloscope
- Connect its tip to the node 1 in the example circuit
- Push *AUTOSET* button
- Push *CH 1 MENU*
- Push the *Coupling* button on the right side of the display and see the various coupling options
- Push *Invert* button to invert polarity of the signal
- Push *CH 2 MENU*

You can also get a plot of one signal versus the other (its called XY mode). Follow the steps below to get a plot with channel-1 voltage on the X-axis and channel-2 voltage on the Y-axis

- Push *CH 1 MENU* and use *Coupling* button to change the coupling to *Ground*
- Use the *POSITION* knob of channel-1 to bring the ground line of channel-1 to the center of the display
- Change the coupling back to *AC*
- Push *CH 2 MENU* and use *Coupling* button to change the coupling to *Ground*
- Use the *POSITION* knob of channel-2 to bring the ground line of channel-2 to the center of the display
- Change the coupling back to *AC*
- Push the *DISPLAY* button
- Push the *FORMAT* button and select *XY* format
- You will see a ellipse shaped plot on the display
- Change the amplitude and frequency of the input sine wave. What do you observe?

Mathematical Operation

DSO has the capability to perform some basic mathematical operations with the input signals.

- Push *MATH MENU* button
- Push the *Operation* button to change the operation to $+$, $-$ and *FFT*